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Bringing into service / operation	Page 33
--	---------

Troubleshooting	Page 34
------------------------------	---------

Cleaning and care	Page 34
--------------------------------	---------

Storage	Page 34
----------------------	---------

Information

Warranty / Service centre.....	Page 34
--------------------------------	---------

Disposal	Page 35
-----------------------	---------

Declaration of conformity	Page 35
--	---------

The following pictograms are used in this operating manual / on the device:

	Please read the operating instructions!		Risk of loss of life by electric shock from damaged mains lead or mains plug!
V~	Volt (AC)		Do not pierce the blanket with needles.
Hz	Hertz (frequency)		Risk of fire!
W	Watt (effective power)		Risk of burns! Warning – Hot surfaces.
	Protection class II		Caution – Risk of electric shock! Risk to life!
	Observe all warnings and safety advice!		The cover 2 is machine-washable. Set the washing machine on a gentle 30 °C washing program.
	Keep children away from the device!		Dispose of the packaging and the device in an environmentally friendly manner!

Shiatsu Chair Cushion

● Introduction

Please make yourself familiar with the device before preparing it for use or using it for the first time. In addition, please refer carefully to the operating instructions below and the important safety advice. Use the device only as described and for the indicated purposes. Keep these instructions in a safe place. Make sure that all documents accompany the device if you pass it on to anyone else.

● Proper use

This Shiatsu cushion is designed for massaging the back. The device is only for private use; it is not intended for medical or commercial purposes and is no substitute for medical treatment. In particular, it must not be used to massage babies, infants,

helpless persons, or animals. The device is only suitable for use in dry indoor locations. Any use other than that described here or modification to the device is not permitted and could lead to injury and / or damage to the product. The manufacturer accepts no liability for any damage caused by improper use. The device is not intended for commercial use.

● Description of parts

Shiatsu cushion:

- 1** Neck cushion
- 2** Cover
- 3** 4 Massage heads (of which 2 are heated and illuminated), rotating vertically in pairs

Hand control:

- 4** ON / OFF button
- 5** Heat button (lighting and heating function)
- 6** A+B button (your whole back)
- 7** A button (upper back)

IMPORTANT ADVICE – KEEP SAFE FOR LATER REFERENCE!

- 8 B button (lower back)
- 9 Timer button
- 10 Button for seat vibration
- 11 Demo button

● Supply Scope

- 1 Shiatsu Chair Cushion
- 1 Cover
- 1 Neck cushion
- 1 Allen key
- 1 Operating instructions

● Technical data

Type:	SSMK 40 B2
Power supply:	220–240 V~ 50/60 Hz, 12 W
Dimensions:	approx. 45 x 104 x 11 cm (W x H x D)
Protection class:	II/□
Electrical field strength:	max. 5000 V/m
Magnetic field strength:	max. 80 A/m
Magnetic flux density:	max. 0.1 Millitesla

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Important safety advice

Failure to comply with the following instructions can lead to personal injury or material damage (electric shocks, skin burns, fire). The following advice on safety and hazards is intended to protect your health and that of other people as well as to protect the device. Please observe these safety instructions and make sure that you pass on this manual together with the device.

⚠ WARNING! DANGER OF LOSS OF LIFE AND ACCIDENT TO INFANTS AND CHILDREN!

Never leave children unsupervised with the packaging materials. Danger of suffocation.

⚠ WARNING! RISK OF ELECTRIC SHOCK!

Ensure that the massage device, switch, mains adapter and lead do not come into contact with water, steam or other liquids. Use the device only in dry indoor environments (e.g. never in the bath, sauna). If a device falls into water, do not attempt to pick it up. Pull the mains adapter out of the mains socket immediately.

Do not use the mains lead for any purpose for which it was not intended, e.g. to carry the product, to hang it up or to pull the mains plug out of the mains socket. Do not allow the lead to become trapped, kinked or twisted. Keep the device and lead away from heat, oil and sharp edges. If the mains cable is damaged or severed, do not touch it, and disconnect the mains plug immediately. Damaged or tangled leads increase the risk of electric shock.

RISK OF ELECTRIC SHOCKS!

Do not insert needles or pointed objects into the device.

Important safety advice

- Do not operate the device during a thunderstorm.
- Operate the device only with the supplied mains adapter and at the mains voltage indicated on the mains adapter.

⚠ WARNING! RISK OF ELECTRIC SHOCK AND INJURY! Do not subject the device to shock and do not drop it. Before every use, check the device carefully for signs of wear and / or damage. Do not use it if you see signs of wear and tear or damage or if the machine has been used incorrectly and switch it off immediately in case of defects or malfunction. In this case, return it to the dealer before using it again.

-  **FIRE HAZARD!** Do not cover the shiatsu massage device during use (cushions, covers etc.). Never use the device near petrol or other readily flammable substances.

-  Ask your doctor and, if relevant, the manufacturer of your heart pacemaker before you use this shiatsu massage device. Under certain circumstances, the electrical and magnetic fields emitted from the device could interfere with the proper functioning of your heart pacemaker. However, the parameters of the fields do not exceed the permitted limits: for information on these limits please refer to the section on "Technical data".

⚠ WARNING! Do not use the device for massaging around the heart if you have a heart pacemaker.

- Children or persons who lack the knowledge or experience to use the device or whose physical, sensory or intellectual capacities are limited must never be allowed to use the device without supervision or instruction by a person responsible for their safety. Children must never be allowed to play with the device.

-  Do not use this device on helpless people, on children, or people who are insensitive to heat (e.g. diabetics, people whose skin is affected by illness or who have scarred skin in the area where the heat will be applied) or animals. Do not use the device during pregnancy, after taking painkilling medication or alcohol.

-  Ensure that you do not fall asleep while the device is operating.

- Do not use the device
 - if you are suffering from a slipped disc or have open wounds,
 - or in cars,
 - not during activities in which an unforeseen reaction could be hazardous.

-  **RISK OF BURNS!** The surface of the device becomes hot during operation. People who are insensitive to heat must be careful while using the device.

- Do not use the shiatsu cushion device on swollen, burned, inflamed, diseased or injured parts or regions of the body or areas of skin. If you are in doubt, seek medical advice before using the device, in particular:

- if you suffer from a serious illness or have had an operation on the upper part of your body,
- if you have a heart pacemaker, implant or other medical device,
- if you are suffering from thrombosis,
- if you suffer from pains with an undiagnosed cause.

-  **CAUTION! DANGER OF CRUSHING!** Ensure that your fingers or other body parts do not get themselves between the rotating massage heads .

- The Shiatsu cushion is designed only massage the backs of persons. It is not a substitute for medical treatment by your doctor.

- The shiatsu cushion device must not be used on the face (e.g. eyes), on the larynx, the feet or other sensitive body parts.

- Repairs may only carried out only by qualified specialist personnel or in a workshop approved by the manufacturer. Unsatisfactory repairs can lead to considerable danger for the user.

⚠ WARNING! If the mains lead to the product is damaged, in order to avoid causing further danger it must be replaced by the manufacturer or his customer services department.

Before use / Bringing into service / operation

● Before use

-  **CAUTION!** Remove the transportation lock for the massage mechanism before you start to use the device.
To do this, remove the screw using the provided tool on the rear side of the cushion. If this screw of the transportation lock is not removed and the Shiatsu cushion is used at one, the cushion will be damaged irreparably and unusable. In this case no warranty will be assumed.
- Always employ the transportation lock screw for the duration of any further transport in the future. Every transport without a transportation lock can damage the device irreparably.

● Bringing into service / operation

Note: Only use the device in an upright position, never lay it flat on the floor.

Note: The zip fastener on the Shiatsu cushion may not be opened because it is only there for technical reasons.

Note: The device may give off an odour of plastic on first use. This odour will cease to develop after a short while.

Note: The device automatically switches off after 15 minutes. Use the shiatsu cushion for a maximum period of 15 minutes. Longer use may lead to over-stimulation and hence tension in the muscles.

Note: Ensure that you are sitting in the middle of the device and that the massage heads [3] can move to the left and right of your spine.

Note: Do not lie, sit and stand with your full weight on the moving parts of the device and do not place any objects on the device.

- Insert the mains plug into the mains socket. Position the lead to avoid creating a trip hazard.

Switching on / off:

- Switch on the Shiatsu cushion by pressing the On / Off button [4].
- After use, switch off the Shiatsu cushion by pressing the On / Off button [4].

Back massage:

- Select with button A+B [6], A [7] or B [8] the area you would like to have massaged:

A+B: your whole back

A: upper back

B: lower back.

The LED lamp signalises the massage area that has been set.

Note: When you select the massage area, the automatic switch-off will be activated and set at 15 minutes.

- Press the ently [9], employing the appropriate frequency, to set the automatic switch-off at 5, 10 or 15 minutes.

Note: By resetting the timer setting, the timer will begin counting anew. At the beginning of the massage, select the maximum duration of the massage and do not change this during the massage in order to prevent the device from overheating.

Note: Ensure that you are experiencing the massage as pleasant and relaxing. If necessary, vary the pressure and/or your position. Stop the massage at once if it is painful or unpleasant.

Heating function:

- Press the heat button [5], to activate the heating function.
- Press the heat button [5] once again to deactivate the heating function.

Note: The heating function cannot be activated unless a massage program has been previously selected.

Seat vibration:

- Press the seat vibration button [10] to activate the seat vibration. The vibration starts at a low intensity and the relevant LED for seat vibration will light up.
- Press the seat vibration button [10], with the appropriate frequency, to switch to low, medium or high intensity or to deactivate the seat vibration.

Note: The seat vibration can be activated together with the back massage. However, the vibration intensity will be slightly reduced in this case.

Demonstration:

- Press the Demo button **[11]** to have the functions of the device demonstrated to you.

● Troubleshooting

- Problem
- Cause
- ◎ Rectification

- Massage heads rotate slower and slower.
- Too much pressure on massage heads.
- ◎ Reduce pressure on massage heads.
- Massage heads do not move.
- Device is not connected to the mains supply.
- ◎ Insert the plug and switch on the device.
- Device is not switched on.
- ◎ Switch on the device by pressing the On/Off button **[4]**. Select the massage area.
- The transportation lock has not been removed. The massage mechanism is damaged.
- ◎ Contact your specialist dealer or the service address.

● Cleaning and care

-  **⚠ WARNING! RISK TO LIFE FROM ELECTRIC SHOCK!** Always disconnect plug from the mains socket before cleaning the shiatsu seat cushion.
-  **RISK OF ELECTRIC SHOCK!** Do not pull the cover **[2]** over the cushion until it is completely dry. There is otherwise a risk of electric shock.
- **CAUTION!** Do not use chemical cleaners or scouring agents for cleaning the shiatsu cushion. Otherwise you may damage it.
- Clean the device only in the manner indicated. Under no circumstances must any liquid be allowed to enter the device or the accessories.
- Small stains can be removed with a cloth or moist sponge and if necessary with a little mild detergent. Do not use solvent-containing cleaners.

-  Clean the removable cover **[2]** in accordance with the cleaning symbols on the label. It is machine-washable. Set the washing machine on a gentle 30 °C washing program.
- Use a mild detergent according to the manufacturer's instructions.

● Storage

- **CAUTION!** Let the shiatsu cushion cool before putting it away.
- If you do not intend to use the Shiatsu cushion for any length of time, store it in the original packing in a dry place without weighing it down.

● Information

● Warranty / Service centre

We provide a 3-year warranty covering faults in materials or manufacture of the product from the purchase date.

The warranty does not cover:

- Damage resulting from improper use.
- Wear parts.
- Defects that were known to the customer at the time of purchase.
- Damage for which the customer was responsible.

This warranty does not affect the customer's legal rights.

The customer must provide proof of purchase in order for any warranty claims within the warranty period to be honoured. Claims under the warranty within a period of 3 years from the date of purchase are honoured by

MGG Elektro GmbH-Service
c/o Dachser GmbH & Co. KG
Thomas-Dachser-Str. 1
89129 Langenau, Germany

Information / Disposal / Declaration of conformity

In the event of a warranty claim, the customer has the right to have the goods repaired by our own workshop or a workshop authorised by use. The device must not be opened for any reason. Opening or modifying the device invalidates the warranty. This warranty does not grant any additional rights to the customer. Complaints are often initiated by difficulties encountered during use. Many of these complaints can be resolved by telephone or e-mail. Please contact our Service Hotline in the first instance, before you return the device to the manufacturer.

GB

☎ 02030 249 050

MT

☎ 800 62036

● Declaration of conformity **CE**

This appliance complies with the requirements of the relevant European and national guidelines. This is confirmed by the CE mark. The relevant declarations are held by the manufacturer.

● Disposal

The packaging is made entirely of recyclable materials, which you may dispose of at local recycling facilities.

Do not dispose of electrical devices with the household rubbish!

In accordance with European Directive 2002/96/EC (covering waste electrical and electronic equipment) and its transposition into national legislation, worn out electrical devices must be collected separately and taken for environmentally compatible recycling.

Contact your local refuse disposal authority for more details of how to dispose of your worn out electrical devices.

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c/o Dachser GmbH & Co.KG